


# Odyssey 1500 Mobile Sink

## PRODUCT SPECIFICATION



tel. 0330 088 3628  
email. [info@mobilesink.com](mailto:info@mobilesink.com)  
[www.mobilesink.com](http://www.mobilesink.com)

### Aesthetics

Unit Base Colour	Brushed Stainless-Steel Grade 304 L
Bowl Colour	Brushed Stainless-Steel Grade 304 L
Control Button Material	Plastic, Standard Black, Red & Blue for hot and cold option

### Function

Adjustable Temp Range	22 to 48 degrees
Adjustable Water Dispense	170ml or 270ml
Water Container Capacity	10 Litres
Removable Water Containers	Yes

### Electrical

Electricity supply Needed	Yes
Voltage (V)	220 – 240 Volts
Frequency Hz	50-60Hz
KW	3KW as standard 2KW on request (2.640 watts and 1.640 watts)
Power Cord	1 Metre, Moulded 3 Pin Plug, 13 amp Fuse

### Dimensions

Footprint (WxD)	320mm x 320mm
Height	850mm
Bowl Depth	150mm
Height to Tap	1040mm
Height, Depth, Length Packed & Palletised	1200mm x 450mm x 450mm
Weight Boxed	18KG

### Technical

Heating System	Thermocoil Technology
Noise Level	55. – 58.0dB (a human voice is around 60dB)
Pump Pressure	2-2.5 Bar
Guarantee	1 Year